COVID-19 is an emerging, rapidly evolving situation. Get the latest public health information from CDC: https://www.coronavirus.gov. Get the latest research from NIH: https://www.nih.gov/coronavirus Find NCBI SARS-CoV-2 literature, sequence, and clinical content: https://www.ncbi.nlm.nih.gov/sars-cov-2/.

Journal List > Elsevier Public Health Emergency Collection > PMC7275188

US National Library of Medicine

National Institutes of Health

Emergency Collection

Public Health

ewlar Public Health Emergency Collection

Esawiar Public Health Emergency Collection

Elseviar Public Health Emergency

Finalgancy

Emargency Collection

Public Health Emergency

Public Health

Emargency Collect

Eseviar Public Health Emergency Collection

Public

Emargency Collection

References

Results

Abstract

Elsevier Public Health Emergency Collection

Public Health Emergency COVID-19 Initiative Am J Infect Control. 2020 Aug; 48(8): 951-954 Published online 2020 Jun 6. doi: 10.1016/j.ajic.2020.06.002

PMCID: PMC7275188 PMID: <u>32522608</u>

Evaluation of an electrostatic spray disinfectant technology for rapid decontamination of portable equipment and large open areas in the era of SARS-CoV-2

Jennifer L. Cadnum, BS,^a Annette L. Jencson, CIC,^a Scott H. Livingston, MD,^b Daniel F. Li, BS,^a Sarah N. Redmond, BS, b Basya Pearlmutter, BS, Brigid M. Wilson, PhD, and Curtis J. Donskey, MDb, MDb, C, *

► Author information ► Copyright and License information <u>Disclaimer</u>

In the setting of the coronavirus disease 2019 pandemic, efficient methods are needed to decontaminate

shared portable devices and large open areas such as waiting rooms. We found that wheelchairs, portable equipment, and waiting room chairs were frequently contaminated with potential pathogens. After minimal manual precleaning of areas with visible soiling, application of a dilute sodium hypochlorite disinfectant using an electrostatic sprayer provided rapid and effective decontamination and eliminated the benign virus bacteriophage MS2 from inoculated surfaces. Key Words: Clostridioides difficile, Wheelchairs, Viruses, Surfaces

Contaminated surfaces are a potential source for dissemination of many bacterial and fungal pathogens. \(^{\dagger} There is also increasing concern that the environment may be an underappreciated source for spread of respiratory viruses, including severe acute respiratory syndrome coronavirus 2.2, 3, 4 Many respiratory viruses survive for hours to days on surfaces, ², ³, ⁴ and respiratory virus nucleic acid has been recovered from surfaces in health care and community settings, including households, day care centers, airports, and

schools. 4, 5, 6, 7 Enhanced environmental cleaning and disinfection is therefore recommended as one component of control measures for severe acute respiratory syndrome coronavirus 2 in health care and community settings.⁸ Most cleaning and disinfection of surfaces is performed through manual application of liquid disinfectants. Cleaning prior to or concordant with application of disinfectants is generally considered important to

reduce bioburden that might decrease the action of the disinfectant. However, in one study, organic and

inorganic material recovered from hospital surfaces did not affect efficacy of sodium hypochlorite and only modestly affected efficacy of ultraviolet-C (UV-C) light. Moreover, thoroughness of manual cleaning is often inadequate, and application can be challenging and time-consuming particularly when surfaces are irregular or large open areas are being cleaned. One potential strategy to improve cleaning and disinfection under these conditions might be to apply disinfectant as a spray with only minimal precleaning to remove visible soil. However, relatively little information is available on this approach to cleaning and disinfection. Here, we tested the effectiveness of a novel spray disinfectant technology that uses an electrostatic sprayer to apply a sporicidal disinfectant to surfaces after minimal precleaning. We focused on items or areas that are challenging to manually clean and disinfect, including wheelchairs, portable equipment, and patient waiting areas. Methods Go to: ₩

Setting The Louis Stokes Cleveland Veterans Affairs Medical Center is a 215-bed acute care facility. Shared-use

wheelchairs are present throughout the facility, and hospital transport staff are available to assist patients.

The facility policy states that wheelchairs should be wiped daily with a disinfectant wipe and as needed if

they become soiled and portable equipment is to be wiped by personnel after each use. Waiting areas are to

be cleaned daily by environmental services personnel, but cleaning of these areas is not monitored.

Test organisms

Test organisms included the nonenveloped single-stranded bacteriophage MS2 (American Type Culture Collection [ATCC] 15597-B1) and Clostridioides difficile spores (ATCC strain 43598). Bacteriophage MS2 was propagated in *Escherichia coli*. ¹⁰ C difficile spores were prepared as previously described. ⁹ Electrostatic sprayer device and disinfectant

Figure 1 shows the electrostatic sprayer device (Clorox Total 360 System-Electrostatic Sprayer, Clorox,

Oakland, CA). The device is 1 m tall and is moved from room to room on wheels. It is plugged into a

standard electrical outlet and a hand-held nozzle is used to direct a fine mist onto surfaces. The sprayer delivers electrostatically charged droplets with an average size of 40-80 µm that are actively attracted to surfaces to improve thoroughness of surface coverage. The device is intended to be used with a variety of different disinfectants and sanitizers. For the current study, the disinfectant used was Spore Defense Cleaner Disinfectant (Clorox) which contains 0.25% sodium hypochlorite. Per the manufacturer, this relatively dilute sodium hypochlorite solution leaves only minimal residual when sprayed and is noncorrosive to common materials in health care settings. No protective equipment is required, but the manufacturer recommends that users wear goggles.



Because the spray bleach product tested has a reduced concentration of sodium hypochlorite in comparison

Go to: ₩

Based on the manufacturer's recommendations, a 5-minute exposure time was used for *C difficile* spores and a 2-minute exposure time was used for MS2; the product has a 1-minute claim against vegetative bacterial pathogens and many viruses. The product was sprayed once at 6 inches from stainless steel carriers applying sufficient disinfectant to thoroughly wet the inoculated surface of the carrier for the

specified exposure time. The carriers were neutralized with Dey-Engley neutralizer (Remel Products, Lenexa, KS). Serial dilutions were plated on selective media and cultured as previously described. 9 Log₁₀ reductions were calculated by subtracting viable organisms recovered after exposure to the disinfectants versus deionized water controls. Experiments were performed in triplicate. Effectiveness in reducing bacterial pathogens and inoculated bacteriophages on surfaces We examined the effectiveness of the spray application of disinfectant on wheelchairs (N=30), portable medical equipment (N=40 devices), and patient waiting area chairs (N=30). The wheelchairs and chairs in waiting areas included soft and hard surfaces. Portable equipment included bladder scanners, electrocardiogram machines, pulse oximeters, workstations on wheels, and Doppler ultrasounds. A commercial improved hydrogen peroxide wipe was used to clean and disinfect areas with obvious visible soiling. The spray disinfectant was then applied once to the surfaces in accordance with the manufacturer's

to many other bleach products, we initially tested its efficacy for killing of C difficile spores (ATCC strain

Agricultural Chemists International Official Method 961.02 Germicidal Spray Products as Disinfectants. 11

43598) and bacteriophage MS2 in the presence of 5% fetal calf serum using Association of Official

longer. Horizontal surfaces typically remained wet for at least 5 minutes but areas at the periphery began to appear visually dry within 2 minutes, and surfaces that were curved or vertically oriented also became dry within 2 minutes. Cultures were collected before spraying and after spraying and air drying. CultureSwabs (Becton Dickinson) premoistened with Dey-Engley neutralizer were used to collect samples from the wheelchairs (composite of arm rest, seat, seat back rest, and hand grips for pushing the wheelchairs), portable equipment (composite including the controls and handles commonly touched during use), and waiting area chairs (composite of chair arm rests, seat, and seat back rest); 10 cm² areas were sampled using a template or one-half of the entire surface of small items such as hand grips were sampled. For C difficile cultures, sterile gloves were donned and sterile 2×2 cm gauze pads premoistened with Dey-Engley neutralizer were used to sample the same sites. After spraying, adjacent 10 cm² areas were sampled or alternate halves of small items were sampled. Cultures were processed for C difficile, methicillin-resistant Staphylococcus

For the wheelchair evaluation, additional work was performed to assess ability to eradicate inoculated

wheelchairs, 10⁶ plaque-forming units (PFU) of bacteriophage MS2 was inoculated and spread to cover 1

disinfectant application. The personnel spraying the disinfectant was blinded to the location where the MS2

viruses and to assess the time required for the spray application versus manual cleaning. For 3 of the

cm² areas of the seat, seat back rest, arm rest, and hand grip and allowed to air dry before the spray

disinfectant was applied; 3 control wheelchairs were inoculated concurrently but not treated with

aureus, enterococci, and gram-negative bacilli as previously described. ⁹

instructions and allowed to air dry. After spraying, all surfaces remained visibly wet for 2 minutes or

was applied. Cultures of swabs and gauze pads for bacteriophage MS2 were processed as previously described. ¹⁰ The time required to apply the spray disinfectant was measured in comparison to the time required for research personnel to manually apply the same disinfectant to thoroughly cover the surfaces on the body of the wheelchair. Data analysis Fisher exact test was used to compare the percentages of contamination of the wheelchairs, portable equipment, and waiting area chairs with a composite of any of the pathogens cultured. All analyses were performed using R version 3.5.1 statistical software (The R Foundation for Statistical Computing, Vienna, Austria).

On steel disks, the sodium hypochlorite spray reduced C difficile spores by $\geq 6.0 \log_{10}$ colony-forming

units with a 5-minute contact time and bacteriophage MS2 by ≥6.0 log₁₀ PFU with a 2-minute contact

time. As shown in Figure 2, contamination with 1 or more of the potential pathogens was present on 30%

or more of surfaces cultured in waiting rooms and on portable equipment and wheelchairs. However, areas

of visible soiling requiring precleaning were uncommon (3 of 30, 10% wheelchairs; 2 of 40, 5% portable

application of the spray disinfectant ($P \le .01$ for each comparison). C difficile contamination was reduced

but not eliminated for each of the sites, whereas the other bacteria were eliminated except for 1 isolate of

devices; 0 of 30 waiting room chairs). There was a significant reduction in contamination for each site after

gram-negative bacilli recovered after spraying a wheelchair. Approximately 5 minutes were required to complete spraying of a waiting room area with 15-20 chairs. 35 30 25 20 20 15 15


■ Before spraying
■ After spraying

(a)

10

■ Before spraying ■ After spraying

(b)



There is no requirement that protective equipment be worn, but the manufacturer does recommend use of goggles during operation. The product left minimal to no residue, whereas disinfectants with higher sodium hypochlorite concentrations often leave a residue. The disinfectant has a 1-minute contact time for killing of many vegetative bacteria and respiratory viruses. One limitation of the dilute sodium hypochlorite product used in this study is that a 5-minute contact time is required for C difficile spores. On real-world surfaces, C difficile spores were reduced but not eliminated completely after a single spray application of the disinfectant. It is likely that failure to eliminate C difficile

approximately 2 minutes. Thus, in settings where C difficile is a concern, repeated application of the dilute

bacteriophage MS2 rather than against viral pathogens. However, there is evidence that bacteriophage MS2

Our study has some limitations. First, we assessed efficacy of the spray disinfectant against the benign

may have increased resistance to liquid disinfectants in comparison to enveloped respiratory viruses. 12

spores was in part related to the fact that curved or vertical surfaces typically had drying times of

sodium hypochlorite product may be required to maintain 5 minutes of wet contact time.

sodium hypochlorite disinfectant using an electrostatic sprayer provided a rapid and effective means to

As noted previously, the electrostatic sprayer can be used with a variety of different disinfectant and

sanitizers. The dilute sodium hypochlorite solution applied in the current study has several advantages.

reduce bacterial contamination on these surfaces and to eliminate an inoculated bacteriophage.

Second, we only studied one type of disinfectant with the spray technology. Third, we did not compare the efficacy and efficiency of the spray technology with alternatives such as UV-C light. However, UV-C light is not well-suited for irregular devices with multiple angles (eg, wheelchairs) that might result in shadowing or large open areas such as waiting rooms. Fourth, we did not collect information on the number of colonies of the pathogens recovered from the surfaces. Finally, we did not compare the efficacy of the spray technology with manual application of disinfectant. However, we do not anticipate that the spray technology will replace manual cleaning of items and surfaces that can easily be wiped. Rather, the technology will be most useful for items and areas that are not amenable to standard cleaning and disinfection.

In summary, our results suggest that application of a dilute sodium hypochlorite disinfectant using an electrostatic sprayer could provide rapid and effective decontamination of portable equipment and large open areas. Additional studies are needed to evaluate the utility of the spray technology in community settings. For example, the technology could be useful for decontamination of areas such as airport waiting

areas, classrooms, and gyms during situations such as the coronavirus disease pandemic.

Footnotes Go to: 🗹 Financial support: This work was supported by a grant from Clorox and by Merit Review grant (1 I01 BX002944-01A1) from the Department of Veterans Affairs to C.J.D. Conflicts of interest: C.J.D. has received research grants from Clorox, Pfizer, and PDI. All other authors report no conflicts of interest relevant to this article.

Go to: ₩

MERS coronaviruses and influenza virus in healthcare settings: possible role of dry surface contamination. J Hosp Infect. 2016;92:235–250. [PMC free article] [PubMed] [Google Scholar] 3. Doremalen N, Bushmaker T, Morris D. Aerosol and surface stability of SARS-CoV-2 as compared with

SARS-CoV-1. N Engl J Med. 2020;382:1564–1567. [PMC free article] [PubMed] [Google Scholar]

2. Otter JA, Donskey C, Yezli S, Douthwaite S, Goldenberg SD, Weber DJ. Transmission of SARS and

1. Donskey CJ. Decontamination devices in health care facilities: practical issues and emerging

applications. Am J Infect Control. 2019;47S:A23–A28. [PubMed] [Google Scholar]

4. Boone SA, Gerba CP. Significance of fomites in the spread of respiratory and enteric viral disease. Appl Environ Microbiol. 2007;73:1687–1696. [PMC free article] [PubMed] [Google Scholar] 5. Boone SA, Gerba CP. The occurrence of influenza A virus on household and day care center fomites. J

Infect. 2005;51:103–109. [PubMed] [Google Scholar] 6. Fong MW, Leung NHL, Xiao J, et al. Presence of influenza virus on touch-surfaces in kindergartens and primary schools [e-pub ahead of print]. J Infect Dis. 10.1093/infdis/jiaa114. Accessed June 24, 2020. [PubMed] [CrossRef]

7. Memish ZA, Almasri M, Assirri A. Environmental sampling for respiratory pathogens in Jeddah airport during the 2013 Hajj season. Am J Infect Control. 2014;42:1266–1269. [PMC free article] [PubMed] [Google Scholar] 8. Centers for Disease Control and Prevention. Infection control guidance for healthcare professionals

9. Zhang A, Nerandzic MM, Kundrapu S, Donskey CJ. Does organic material on hospital surfaces reduce the effectiveness of hypochlorite and UV radiation for disinfection of Clostridium difficile? Infect Control Hosp Epidemiol. 2013;34:1106–1108. [PubMed] [Google Scholar]

10. Tomas ME, Kundrapu S, Thota P. Contamination of the skin and clothing of healthcare personnel

during removal of personal protective equipment. JAMA Intern Med. 2015;175:1904–1910. [PubMed]

about coronavirus (COVID-19). Available at: https://www.cdc.gov/coronavirus/2019-ncov/index.html.

[Google Scholar] 11. AOAC Official Method 961.02 Germicidal spray products as disinfectants. Microchem Laboratory website. Available at: http://microchemlab.com/test/aoac-germicidal-spray-products-test-aoac-96102.

2015. Accessed April 10, 2020. 12. Sattar SA. Hierarchy of susceptibility of viruses to environmental surface disinfectants: a predictor of activity against new and emerging viral pathogens. J AOAC Int. 2007;90:1655–1658. [PubMed]

Formats:

Go to: 🗹

Links

Article | PubReader | ePub (beta) | PDF (701K) | Citation Share 🗗 Facebook 💟 Twitter 🛭 👺 Google+

Sign in to NCBI

Help

Save items

Add to Favorites Similar articles in PubMed Ultraviolet Germicidal Irradiation to Decontaminate Filtering Face

Piece Respirators During COVID-19 Pandemic. [S D Med. 2020] The ABC's of disaster management: Managing apheresis operations during the SARS-CoV-2 pandemic. [J Clin Apher. 2020] Protection and disinfection policies against SARS-CoV-2 (COVID-[Infez Med. 2020] Effectiveness of N95 Respirator Decontamination and Reuse

against SARS-CoV-2 Virus. [Emerg Infect Dis. 2020] Evidence for decontamination of single-use filtering facepiece respirators. [J Hosp Infect. 2020] See reviews.. See all.

PubMed **Recent Activity** Turn Off Clear

Evaluation of an electrostatic spray disinfectant technology for rapid decontami... ECG Diagnosis: Complete Heart Block **Biocidal Effectiveness of Selected Disinfectants Solutions**

Identification of Escherichia coli from broiler chickens in Jordan, their antimi... See more..

Based on Water and Oz...

[Google Scholar]

Accessed April 10, 2020.